

LOGIQ S7 Expert

Amazing **versatility**

Data Sheet



➤ **Sensational** performance

➤ **Smart** design

➤ **Specialized** capabilities

➤ **Product Description**

The LOGIQ* S7 Expert is a highly mobile and easy to use, performance multipurpose color Doppler imaging system, designed for Obstetrics, Gynecology, Cardiology, Vascular, Urology, Small Parts, Pediatric, Neonatal, Transcranial, and Abdominal applications.



General specifications

Dimensions and weight

Height	Standard: 1750 mm (68.9 in) Tall: 1115 mm (43.9 in)
Width	Keyboard: 500 mm (19.7 in) Caster: 620 mm (24.4 in)
Depth	Maximum: 856 mm (33.7 in) Caster: 790 mm (31.1 in)
Weight (no Peripherals)	90 kg/198 lbs

Electrical power

Voltage	100-120 Vac or 220-240 Vac
Frequency	50/60 Hz
Power consumption maximum of 900 VA with peripherals	

Console design

4 active probe ports, 1 non-imaging
Integrated HDD and DVD-R/W
On-board storage for peripherals
Integrated speakers
Probe holders
Gel holder/warmer
Front and rear handles

User interface

Operator keyboard

Ergonomic full size keyboard
Swivel-adjustable, Height-adjustable
8 TGC pods
7" (177.8 mm) wide LCD touch screen

Monitor

19" (482.6 mm) high-resolution LCD
Articulating monitor arm

System overview

Applications

Abdominal	Obstetrical
Gynecological	Breast
Small Parts and Superficial	Musculoskeletal
Vascular	Urological
Endocavitary	Pediatric and Neonatal
Transcranial	Cardiac

Scanning methods

Electronic Sector
Electronic Convex
Electronic Micro Convex
Electronic Linear
Real Time 4D Volume Sweep

Transducer types

Sector Phased Array
Convex Array
Microconvex Array
Linear Array
Matrix Array
Single CW (Pencil) Probes
Volume Probes (4D)

Operating modes

B-Mode
Coded Harmonic Imaging
M-Mode
Color Flow Mode (CFM)
Power Doppler Imaging (PDI)
PW Doppler with High PRF
M-Color Flow Mode
Anatomical M-Mode
Curved Anatomical M-Mode
B-Flow/B-Flow Color (Option)
Extended Field of View (LOGIQView Option)
Coded Contrast Imaging (Option)
CW Doppler Mode (Option)
TVI Mode (Option)
Elastography (Option)
3D/4D Volume Modes (Option)

System standard features

Advanced user interface with high resolution 7" wide LCD touch screen
Automatic Optimization
CrossXBeam compounding
Speckle Reduction Imaging (SRI-HD)
Fine Angle Steering
Coded Harmonic Imaging
Virtual Convex
Patient information Database
Image Archive on integrated CD/DVD and hard drive

Raw Data Analysis
Real-time automatic Doppler calculations
OB Calculations
Fetal Trending
Multigestational Calculations
Hip Dysplasia Calculations
Gynecological Calculations
Vascular Calculations
Urological Calculations
Renal Calculations
Cardiac Calculations
Remote capability: InSite ExC
On-board electronic documentation
MPEGVue
Key Macro
Network Storage
Quick Save
Quick Patient Entry

System Options

Auto IMT
Elastography
Elastography Q-Analysis ¹
Advanced 3D
DICOM 3.0 Connectivity
LOGIQView
B-Flow/B-Flow Color
CF/PDI Quantification
B Steer+
Stress Echo
Tissue Velocity Imaging (TVI) with Q-Analysis
Scan Assistant
Report Writer
Coded Contrast Imaging ²
ECG + AHA/IEC Cables
CW Doppler
DVR Kit
Real Time 4D
4D TUI
VOCAL
VCI Static

Cabinet: High/Mid/Low
Drawer
Small Probe Adaptor
Vertical Endocavitary Probe Holder
Side Probe Holder
Probe Cable Hanger
3-Pedal Foot Switch
Isolation transformer

Peripheral Options

Integrated Options for
• Digital BW thermal printer
• Digital A5 Color thermal printer
• DVD video recorder
Digital A6 Color thermal printer
External USB printer connection
HDMI output available for compatible devices
Foot Switch with programmable functionality
Console Protective Cover

Display modes

Live and Stored Display Format: full size and split screen – both with thumbnails for still and Cine
Review Image Format: 4x4 and “thumbnails” for still and Cine
Simultaneous Capability
B or CrossXBeam/PW
B or CrossXBeam/CFM or PDI
B/M
B/CrossXBeam
Real-time Triplex Mode (B or CrossXBeam + CFM or PDI/PW or CW (Option))
Selectable alternating Modes
B or CrossXBeam/PW
B or CrossXBeam + CFM (PDI)/PW(CW(Option))
B/CW (Option)
Multi-image (split/quad screen)
Live and/or frozen
B or CrossXBeam + B or CrossXBeam/CFM or PDI
Independent Cine playback
Time line display
Independent Dual B or CrossXBeam/PW Display

CW
Display Formats
<div>• Top/Bottom selectable format</div>
<div>• Side/Side selectable format</div>
Virtual Convex
Timeline only
Display annotation
Patient Name: First, Last and Middle
Patient ID
2 nd Patient ID
Age, Sex and Birth Date
Hospital Name
Date format: 3 Types selectable
<div>• MM/DD/YY</div> <div>• DD/MM/YY</div> <div>• YY/MM/DD</div>
Time format: 2 types selectable
<div>• 24 hours</div> <div>• 12 hours</div>
Gestational Age from
<div><div>• LMP</div><div>• GA</div></div> <div><div>• EDD</div><div>• BBT</div></div>
Displayed Acoustic Output
<div>• TIS: Thermal Index Soft Tissue</div> <div>• TIC: Thermal Index Cranial (Bone)</div> <div>• TIB: Thermal Index Bone</div> <div>• MI: Mechanical Index</div>
% of Maximum Power output
Probe Name
Map Names
Probe Orientation
Depth Scale Marker
Lateral Scale Marker
Focal Zone Markers
Image Depth
Zoom Depth
B-Mode
Gain
Dynamic Range
Imaging Frequency
Frame Averaging
Acoustic Frame Rate
Gray Map
SRI-HD
M-Mode

Gain
Dynamic Range
Time Scale
Doppler Mode
Gain
Angle
Sample Volume Depth and Width
Wall Filter
Velocity and/or Frequency Scale
Spectrum Inversion
Time Scale
PRF
Doppler Frequency
Color Flow Mode
Line Density
Frame Averaging
Packet Size
Color Scale: 3 types
<div>• Power</div> <div>• Directional PDI</div> <div>• Symmetrical Velocity Imaging</div>
Color Velocity Range and Baseline
Color Threshold Marker
Color Gain
PDI
Inversion
Doppler Frequency
TGC Curve
Cine Gage, Image Number/Frame Number
Body Pattern: Multiple human and animal types
Application Name
Measurement Results
Operator Message
Biopsy Guide Line and Zone
Heart Rate

General System Parameters

System Setup
Pre-programmable Categories
User Programmable Preset Capability
Factory Default Preset Data
Languages: English, French, German, Spanish, Italian, Portuguese, Russian, Greek, Swedish, Danish, Dutch, Finnish,

Norwegian, Japanese (message only)
OB Report Formats including Tokyo Univ., Osaka Univ., USA, Europe, and ASUM
User Defined Annotations
Body Patterns
Customized Comment Home Position

Complete User Manual available on-board through Help (F1)
User Manual and Service Manual are included on CD with each system. A printed manual is available upon request.

CINE Memory/Image Memory
384 MB of Cine Memory
Selectable Cine Sequence for Cine Review
Prospective Cine Mark
Measurements/Calculations and Annotations on Cine Playback
Scrolling timeline memory
Dual Image Cine Display
Quad Image Cine Display
Cine Gauge and Cine Image Number Display
Cine Review Loop
Cine Review Speed

Image Storage
On-board database of patient information from past exams
Storage Formats:
<div>• DICOM – compressed/uncompressed, single/multiframe, with/without Raw Data</div> <div>• Export JPEG, JPEG2000, WMV (MPEG 4) and AVI formats</div>
Storage Devices:
<div>• USB Memory Stick: 64MB to 4GB (for exporting individual images/clips)</div> <div>• CD-RW storage: 700MB</div> <div>• DVD storage: -R (4.7GB)</div> <div>• Hard Drive Image Storage: ~112GB</div>
Compare old images with current exam
Reload of archived data sets

Connectivity & DICOM
Ethernet network connection
DICOM 3.0 (Optional)
Verify
Print
Store
Modality Worklist

Storage Commitment
Modality Performed Procedure Step (MPPS)
Media Exchange
Off network/mobile storage queue
Query/Retrieve
Public SR Template
<div>• Structured Reporting – compatible with vascular and OB standard</div>
Remote capability InSite ExC

Physiological Input Panel (Option)
Physiological Input
ECG, 2 lead
Dual R-Trigger
Pre-settable ECG R Delay Time
Pre-settable ECG Position
Adjustable ECG Gain Control
Automatic Heart Rate Display

Report Writer (Option)
On-board reporting package automates report writing
Formats various exam results into a report suitable for printing or reviewing on a standard PC
Exam result reports can include patient info, exam info, measurements, calculations, images, comments and physician diagnosis
Standard templates provided
Customizable templates

Scanning Parameters
Displayed Imaging Depth: 0 – 33 cm
Minimum Depth of Field: 0 – 2 cm (Zoom) (probe dependent)
Maximum Depth of Field: 0 – 33 cm (probe dependent)
Continuous Dynamic Receive Focus/Continuous Dynamic Receive Aperture
Adjustable Dynamic Range
Adjustable Field of View (FOV)
Image Reverse: Right/Left
Image Rotation of 0°, 180°

Digital B-Mode

Adjustable: <ul style="list-style-type: none">• Acoustic Power• Gain• Dynamic Range• Frame Averaging• Gray Scale Map• Frequency• Line Density• Scanning Size (FOV or Angle – depending on the probe, see probe specifications)• B Colorization• Reject• Suppression• SRI-HD• Edge Enhance

Digital M-Mode

Adjustable: <ul style="list-style-type: none">• Acoustic Power• Gain• Dynamic Range• Gray Scale Map• Frequency• Sweep Speed• M Colorization• M Display Format• Rejection
--

Anatomical M-Mode

M-Mode cursor adjustable at any plane
Can be activated from a Cine loop from a live or stored image
M and A capability
Available with Color Flow Mode
Curved Anatomical M-Mode

Digital Spectral Doppler Mode

Adjustable: <ul style="list-style-type: none">• Acoustic Power• Gain• Dynamic Range• Gray Scale Map• Transmit Frequency• Wall Filter• PW Colorization• Velocity Scale Range• Sweep Speed• Sample Volume Length• Angle Correction• Steered Linear• Spectrum Inversion
--

<ul style="list-style-type: none">• Trace Method• Baseline Shift• Doppler Auto Trace• Time Resolution• Compression• Trace Direction• Trace Sensitivity
--

Digital Color Flow Mode

Adjustable: <ul style="list-style-type: none">• Acoustic Power• Color Maps, including velocity-variance maps• Gain• Velocity Scale Range• Wall Filter• Packet Size• Line Density• Spatial Filter• Steering Angle• Baseline Shift• Frame Average• Threshold• Accumulation mode• Sample Volume Control• Flash Suppression• Quantification (Option)

Digital Power Doppler Imaging

Adjustable: <ul style="list-style-type: none">• Acoustic Power• Color Maps including velocity-variance maps• Gain• Velocity Scale Range• Wall Filter• Packet Size• Line Density• Spatial Filter• Steering Angle• Frame Average• Threshold• Accumulation mode• Sample Volume Control• Flash Suppression

Continuous Wave Doppler (Option)

Adjustable: <ul style="list-style-type: none">• Acoustic Power• Gain• Dynamic Range• Gray Scale Map• Transmit Frequency

<ul style="list-style-type: none">• Wall Filter• CW Colorization• Velocity Scale Range• Sweep Speed• Angle Correction• Spectrum Inversion• Trace Method• Baseline Shift• Doppler Auto Trace• Compression• Trace Direction• Trace Sensitivity

Automatic Optimization

Optimize B-Mode image to improve contrast resolution
Selectable amount of contrast resolution improvement (low, medium, high)
Auto-Spectral Optimize adjusts <ul style="list-style-type: none">• Baseline• Invert• PRF (on live image)• Angle correction

Coded Harmonic Imaging

Available on all 2D probes

B-Flow (Option)

Available on C1-5-D, 9L-D, ML6-15, 11L-D and L8-18i-D probes
Background: On/Off
Sensitivity/PRI
Line Density
Edge Enhance
Frame Average
Gray Scale Map
Tint Map
Dynamic Range
Rejection
Gain
Dual Beam
B-Flow Color
Accumulation

Coded Contrast Imaging (Option)

Available on C1-5-D probe
2 Contrast Timers
Timed Updates: 0.05 – 10 seconds

Accumulation mode, six levels
Maximum Enhance Mode
Flash
Time Intensity Curve (TIC) Analysis
Auto MI control

The LOGIQ S7 Expert is designed for compatibility with commercially available ultrasound contrast agents. Because the availability of these agents is subject to government regulation and approval, product features intended for use with these agents may not be commercially marketed nor made available before the contrast agent is cleared for use. Contrast related product features are enabled only on systems for delivery to an authorized country or region of use.

LOGIQView (Option)

Extended Field of View Imaging
Available on 9L-D, ML6-15, 11L-D, L8-18i-D, 3CRF-D, C1-5-D, IC5-9-D, 3Sp-D, RAB4-8-D, 8C and S4-10-D probes
For use in B-Mode
CrossXBeam is available on linear probes
Auto detection of scan direction
Pre or post-process zoom
Rotation
Auto fit on monitor
Measurements in B-Mode

3D

Allows unlimited rotation and planar translations
3D reconstruction from Cine sweep

Advanced 3D (Option)

Acquisition of Color data
Automatic rendering
3D Landscape technology
3D Movie

Scan Assistant (Option)

Factory Programs
User defined programs
Steps include image annotations, mode transitions, basic imaging controls and measurement initiation

Elastography (Option)

Available on ML6-15, 9L-D, C1-5-D, IC5-9-D, and 11L-D probes
Semi-Quantification ¹

TVI (Option)

Myocardial Doppler Imaging with color overlay on tissue image
Available on the sector probes
Tissue color overlay can be removed to show just the 2D image, still retaining the tissue velocity information
Curved Anatomical M-Mode: free (curved) drawing of M-Mode generated from the cursor independent from the axial plane
Q-Analysis: Multiple Time Motion trace display from selected points in the myocardium

Stress Echo (Option)

Advanced and flexible Stress Echo examination capabilities
Provides exercise and pharmacological protocol templates
8 default templates
Template editor for user configuration of existing templates or creation of new templates
Reference scan display during acquisition for stress level comparison (dual screen)
Baseline level/Previous level selectable
Raw data continuous capture
Over 100 sec available
Wall motion scoring (bulls-eye and segmental)
Smart stress: Automatically set up various scanning parameters (for instance, geometry, frequency, gain etc.) according to same projection on previous level

Virtual Convex

Provides a convex field of view
Compatible with CrossXBeam
Available on all linear and sector transducers

SRI-HD

Speckle Reduction Imaging
Provides multiple levels of speckle reduction
Compatible with Side by Side DualView Display
Compatible with all linear, convex and sector transducers
Compatible with B-Mode, Color, Contrast Agent and 3D imaging

CrossXBeam

Provides 3, 5, 7 or 9 angles of spatial compounding
Live Side by Side DualView Display

Compatible with:

- Color Mode
- PW
- SRI-HD
- Coded Harmonic Imaging
- Virtual Convex

Available on 9L-D, ML6-15, 11L-D, L8-18i-D, 3CRF-D, C1-5-D, RAB4-8-D, 8C and IC5-9-D probes

Controls Available While “Live”

Write Zoom
B/M/CrossXBeam-Mode
Gain
TGC
Dynamic Range
Acoustic Output
Transmission Focus Position
Transmission Focus Number
Line Density Control
Sweep Speed for M-Mode
Number of Angles for CrossXBeam
PW-Mode
Gain
Dynamic Range
Acoustic Output
Transmission Frequency
PRF
Wall Filter
Spectral Averaging
Sample Volume Gate <ul style="list-style-type: none">• Length• Depth
Velocity Scale
Color Flow Mode
CFM Gain
CFM Velocity Range
Acoustic Output
Wall Echo Filter
Packet Size
Frame Rate Control
CFM Spatial Filter

CFM Frame Averaging
CFM Line Resolution
Frequency/Velocity Base Line Shift

Controls Available on “Freeze” or Recall

Automatic Optimization
SRI-HD
CrossXBeam – Display non-compounded and compounded image simultaneously in split screen
3D reconstruction from a stored Cine loop
B/M/CrossXBeam Mode
Gray Map Optimization
TGC
Colorized B and M
Frame Average (loops only)
Dynamic Range: Anatomical M-Mode
Max Read Zoom to 8x: Base Line Shift
Sweep Speed
PW Mode
Gray Map
Post Gain
Baseline shift
Sweep Speed
Invert Spectral wave form
Compression
Rejection
Colorized Spectrum
Display Format
Doppler Audio
Angle Correct
Quick Angle Correct
Auto Angle Correct
Color Flow
Overall Gain (loops and stills)
Color Map
Transparency Map
Frame Averaging (loops only)
Flash Suppression
CFM Display Threshold
Spectral Invert for Color/Doppler
Anatomical M-Mode on Cine loop

Measurements/Calculations

General B-Mode

Depth and Distance
Circumference (Ellipse/Trace)
Area (Ellipse/Trace)
Volume (Ellipsoid)
% Stenosis (Area or Diameter)
Angle between two lines

General M-Mode

M-Depth
Distance
Time
Slope
Heart Rate

General Doppler Measurements/Calculations

Velocity
Time
A/B Ratio (Velocities/Frequency Ratio)
PS (Peak Systole)
ED (End Diastole)
PS/ED (PS/ED Ratio)
ED/PS (ED/PS Ratio)
AT (Acceleration Time)
ACCEL (Acceleration)
TAMAX (Time Averaged Maximum Velocity)
Volume Flow (TAMEAN and Vessel Area)
Heart Rate
PI (Pulsatility Index)
RI (Resistivity Index)

Real-time Doppler Auto Measurements/Calculations

PS (Peak Systole)
ED (End Diastole)
MD (Minimum Diastole)
PI (Pulsatility Index)
RI (Resistivity Index)
AT (Acceleration Time)
ACC (Acceleration)
PS/ED (PS/ED Ratio)
ED/PS (ED/PS Ratio)

HR (Heart Rate)
TAMAX (Time Averaged Maximum Velocity)
PVAL (Peak Velocity Value)
Volume Flow (TAMEAN and Vessel Area)

OB Measurements/Calculations

Gestational Age by: <ul style="list-style-type: none">GS (Gestational Sac)CRL (Crown Rump Length)FL (Femur Length)BPD (Biparietal Diameter)AC (Abdominal Circumference)HC (Head Circumference)APTD x TTD (Anterior/Posterior Trunk Diameter by Transverse Trunk Diameter)FTA (Fetal Trunk Cross-sectional Area)HL (Humerus Length)BD (Binocular Distance)FT (Foot Length)OFD (Occipital Frontal Diameter)TAD (Transverse Abdominal Diameter)TCD (Transverse Cerebellum Diameter)THD (Thorax Transverse Diameter)TIB (Tibia Length)ULNA (Ulna Length)	
Estimated Fetal Weight (EFW) by: <ul style="list-style-type: none">AC, BPDAC, BPD, FL, HCAC, FL, HCBPD, APTD, TTD, FL	
Calculations and Ratios <ul style="list-style-type: none">FL/BPDFL/HCCI (Cephalic Index)CTAR(Cardio-Thoracic Area Ratio)	<ul style="list-style-type: none">AC, BPD, FLAC, FLAC, HCBPD, APTD, TTD, SLFL/ACHC/ACAFI (Amniotic Fluid Index)
Measurements/Calculations by: ASUM, ASUM 2001, Berkowitz, Bertagnoli, Brenner, Campbell, CFEF, Chitty, Eik-Nes, Ericksen, Goldstein, Hadlock, Hansmann, Hellman, Hill, Hohler, Jeanty, JSUM, Kurtz, Mayden, Mercer, Merz, Moore, Nelson, Osaka University, Paris, Rempen, Robinson, Shepard, Shepard/Warsoff, Tokyo University, Tokyo/Shinozuka, Yarkoni	
Fetal Graphical Trending	
Growth Percentiles	
Multi-Gestational Calculations (4)	
Fetal Qualitative Description (Anatomical survey)	
Fetal Environmental Description (Biophysical profile)	
Programmable OB Tables	
Over 20 selectable OB Calculations	
Expanded Worksheets	

GYN Measurements/Calculations

Right Ovary Length, Width, Height
Left Ovary Length, Width, Height
Uterus Length, Width, Height
Cervix Length, Trace
Ovarian Volume
ENDO (Endometrial thickness)
Ovarian RI
Uterine RI
Follicular measurements
Summary Reports

Vascular Measurements/Calculations

SYS DCCA (Systolic Distal Common Carotid Artery)
DIAS DCCA (Diastolic Distal Common Carotid Artery)
SYS MCCA (Systolic Mid Common Carotid Artery)
DIAS MCCA (Diastolic Mid Common Carotid Artery)
SYS PCCA (Systolic Proximal Common Carotid Artery)
DIAS PCCA (Diastolic Proximal Common Carotid Artery)
SYS DICA (Systolic Distal Internal Carotid Artery)
DIAS DICA (Systolic Distal Internal Carotid Artery)
SYS MICA (Systolic Mid Internal Carotid Artery)
DIAS MICA (Diastolic Mid Internal Carotid Artery)
SYS PICA (Systolic Proximal Internal Carotid Artery)
DIAS PICA (Diastolic Proximal Internal Carotid Artery)
SYS DECA (Systolic Distal External Carotid Artery)
DIAS DECA (Diastolic Distal External Carotid Artery)
SYS PECA (Systolic Proximal External Carotid Artery)
DIAS PECA (Diastolic Proximal External Carotid Artery)
VERT (Systolic Vertebral Velocity)
SUBCLAV (Systolic Subclavian Velocity)
Automatic IMT
Summary Reports

Urological Calculations

Bladder Volume
Prostate Volume
Lt/Rt Renal Volume
Generic Volume
Post-Void Bladder Volume

Probes

3CRF-D

Micro Convex Biopsy Probe	
Applications	Abdomen, OB/GYN, Urology
Biopsy Guide	Single-Angle, disposable with a reusable bracket (40442LR), Multi-Angle with a reusable bracket (H40452LP)

8C

Micro Convex Probe	
Applications	Neonatal, Pediatrics
Biopsy Guide available	None

C1-5-D

Convex Probe	
Applications	Abdomen, Vascular, OB/Gyn, Urology
Biopsy Guide	Multi-Angle, disposable with a reusable bracket (H40432LE)

IC5-9-D

Endo Micro Convex Probe	
Applications	OB/GYN, Urology, Transvaginal, Transrectal
Biopsy Guide	Single Angle, disposable with a disposable bracket (E8385MJ, E8333JB), Reusable bracket (H40412LN)

ML6-15

Matrix Array Linear Probe	
Applications	Small parts, Vascular, Neonatal, Pediatrics
Biopsy Guide	Multi-Angle, disposable with a reusable bracket (H40432LJ)

11L-D

Linear Probe	
Applications	Vascular, Small Parts, Neonatal, Pediatrics
Biopsy Guide	Multi-Angle, disposable with a reusable bracket (H40432LC)

9L-D

Linear Probe	
Applications	Vascular, Small Parts, Pediatric, Abdomen
Biopsy Guide	Multi-Angle, disposable with a reusable bracket (H4906BK)

L8-18i-D

Linear Probe	
Applications	Vascular, Small Parts, Neonatal, Pediatrics
Biopsy Guide available	None

3Sp-D

Phased Array Sector Probe	
Applications	Cardiac, Transcranial, Abdomen
Biopsy Guide	Multi-Angle, Reusable bracket (H46222LC)

S4-10-D

Phased Array Sector Probe	
Applications	Pediatrics, Neonatal, Abdomen
Biopsy Guide available	None

RAB4-8-D

Convex Volume Probe	
Applications	Abdomen, OB/GYN, Urology
Biopsy Guide	Single-Angle, disposable with reusable bracket (H46701AE), single angle reusable (H48621Y)

P2D

CW Split Crystal Probe	
Applications	Cardiac, Pediatric

P6D

CW Split Crystal Probe	
Applications	Cardiac, Pediatric, Vascular

Inputs and Outputs

HDMI Out
Ethernet Network (RJ45)
External Audio Out
USB (2x in front, 1x in rear)
AC Power Input

Safety Conformance

The LOGIQ S7 Expert is:

Classified to UL 60601-1 by a Nationally Recognized Test Lab
Certified to CAN/CSA-C22.2 No. 601.1-M90 by an SCC accredited Test Lab
CE Marked to Council Directive 93/42/EEC on Medical Devices
Conforms to the following standards for safety:
<ul style="list-style-type: none">• IEC 60601-1 Medical electrical equipment – Part 1: General requirements for safety• IEC 60601-1-1 Medical electrical equipment – Part 1-1 General requirements for safety – Collateral Standard: Safety requirements for medical electrical systems• IEC 60601-1-2 Medical electrical equipment – Part 1-2 General requirements for safety – Collateral Standard: Electromagnetic compatibility – requirements and tests

- IEC 60601-1-4 Medical electrical equipment Part 1-4 General requirements for safety – Collateral Standard: programmable electrical medical systems
 - IEC 60601-1-6 Medical electrical equipment Part 1-6 General requirements for basic safety and essential performance – Collateral Standard: Usability
 - IEC 60601-2-37 Medical electrical equipment – Part 2-37: Particular requirements for the safety of ultrasonic medical diagnostic and monitoring equipment
 - ISO 10993-1 Biological evaluation of medical devices – Part 1 Evaluation and testing
 - NEMA UD2 Acoustic output measurement standard for diagnostic ultrasound equipment
 - NEMA UD3 Standard for real time display of thermal and mechanical acoustic output indices on diagnostic ultrasound equipment (MI, TIS, TIB, TIC)
- EMC Emissions Group 1 Class B device requirements as per Sub clause 4.2 of CISPR 11

GE Healthcare
9900 Innovation Drive
Wauwatosa, WI 53226
U.S.A.

www.gehealthcare.com

APAC
T +81 42 585 5111

CHINA
T +86 21 3877 4045

EAGM
T +971 4 4296161
or +971 4 4296101

EUROPE
T +49 212 28 020

INDIA
T +91 80 41801000

LATIN AMERICA
T 0 800 122 345

USA/CANADA
T +1 888 202 5528

¹Elastography with semi-Quantification (Elastography Q-Analysis) described in this material has not been cleared by the U.S. FDA and is not available for promotion or sale in the United States.

²Coded contrast imaging described in this material has not been cleared by the U.S. FDA and is not available for promotion or sale in the United States.

©2012 General Electric Company – All rights reserved.

General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Representative for the most current information.

GE and GE Monogram are trademarks of General Electric Company.

GE Medical Systems Ultrasound & Primary Care Diagnostics, LLC, a General Electric company, doing business as GE Healthcare.

* Trademark of General Electric Company.



imagination at work